COVID-19 Ag

INTENDED USE

SGTi-flex COVID-19 Ag is an immunoassay for qualitative detection of Nucleocapsid protein antigen from SARS-CoV-2 in nasopharyngeal and nasal swab specimens. The test is used as an aid in the rapid diagnosis of SARS-CoV-2 viral infections.

SUMMARY AND EXPLANATION

The novel coronavirus (SARS-CoV-2) was identified in December 2019, and in February 2020, the World Health Organization (WHO) officially named the disease caused by SARS-CoV-2 as COVID-19 (Coronavirus Disease 2019). Belonging to the family Coronaviridae, it has a positive-sense single-stranded RNA and can be transmitted between people. The coronaviruses identified for human infection include 229E, NL63 belonging to α-Coronaviruses and HKU1, OC43, SARS-CoV, MERS-CoV belonging to β-Coronaviruses.

The new coronavirus was published under the name of SARS-CoV-2, with 80% of genetic similarity to SARS-CoV by ICTV (International Committee on Taxonomy of Viruses).

COVID-19 spreads mainly through respiratory droplets, which cause lethargy, fever, dry cough, and dyspnea when infected. It can be even led to death with its severe symptoms like sepsis, MOF (Multiple Organ Failure) and ARDS (Acute Respiratory Distress Syndrome). It is more contagious than SARS which caused more than 800 deaths and 8,000 infected patients. Moreover, it has an incubation period of about 3 days to up to 16 days and becomes a big threat as infectivity appears even during the incubation period. There is currently no specific treatment for COVID-19, and rapid and accurate diagnosis is an important issue for isolation of patients with symptoms of suspected COVID-19.

PRINCIPLE

SGTi-flex COVID-19 Ag is an immunoassay for qualitative detection of SARS-CoV-2 antigens directly from nasopharyngeal and nasal swab specimens. The SARS-CoV-2 antigens are extracted from swab in the extraction buffer and the extracted sample solutions are loaded to the sample well of the Test Cassette. When the sample is loaded, the detection antibody binds to SARS-CoV-2 antigen and flows through the membrane. The detection antibody-gold conjugate and SARS-CoV-2 antigen move to the test line area and are accumulated by the capture antibody immobilized on the membrane. This leads to the generation of a reddish colored band. The intensity of the band depends on quantity of SARS-CoV-2 antigen and the test results are interpreted by user's eye according to the instructions for use.

MATERIALS SUPPLIED (CAGT001EO)

- Test Cassette 1
- Extraction Buffer 1 (0.3 mL/tube)
- Dropping cap 1
- Sample collection swab 1
- Instructions for Use 1

STORAGE AND STABILITY

- Store SGTi-flex COVID-19 Ag Test Cassette and Extraction Buffer at 2~30°C (36~86°F).
- If SGTi-flex COVID-19 Ag Test Cassette and Extraction Buffer are stored in cold storage, allow them for 30 minutes to return to room temperature before testing.
- Do not open the pouch of Test Cassette until ready to use. After opening aluminum pouch, Test Cassette should be used immediately.
- Keep away from direct sunlight.

WARNING AND PRECAUTIONS

- For *in vitro* diagnostic use only.
- Clinical diagnosis through this product should be made through a comprehensive review of the specialist based on other test methods and clinical symptoms.
- Please read the instructions carefully before you begin the test and follow the procedure correctly.
- It is prohibited to reuse Test Cassettes because they are single use only.
- The test result after the expiry date is not reliable.
- Test Cassette is sensitive to moisture and should be stored in a sealed pouch until use. Use Test Cassette immediately after opening the pouch.
- Do not use the Test Cassette if it is broken or the pouch is not stored in sealed.
- Samples and Test Cassette must be at room temperature before testing.
- It is an *in vitro* diagnostic product and the risk of infection is low because there is no direct contact with the human body. However please be cautious when handling Test Cassette and samples because of the use of clinical samples containing potential infectious sources. Dispose of the used samples and Test Cassettes properly in accordance with the relevant regulations.
- Smoking and eating are prohibited at test site when handing specimens or kit reagents.

TEST PREPARATION
1. Test should be done immediately after sample collecting.

1) If sample swabs are not used immediately after sample collection, specimen is recommended to be stored in deep freezer at -70°C (or in dry ice or liquid nitrogen). A freezer at -20°C is NOT recommended. If the specimen is stored at 2-8°C, it can be stored up to 72 hours.

2. Preparation before Test

- 1) All samples and reagents should be stored at room temperature and stayed homogenous 15~30 minutes prior to testing.
- 2) Test cassette is moisture sensitive so should be used **immediately** after opening.

SAMPLE COLLECTION

SGTi-flex COVID-19 Ag can be performed with nasopharyngeal swab and nasal swab.


1. Before sample collection, prepare the extraction buffer.
Remove the sealing foil from the Extraction Buffer tube.
2. Place the Extraction Tube in the packing box.


3. SGTi-flex COVID-19 Ag uses the sample of nasopharyngeal swab and nasal swab.

(1) Direct Swab


- 1) Please use single use sample collecting swab.
- 2) For a nasal swab, carefully insert the entire absorption tip of the swab (usually 1/2 to 3/4 of an inch (1 to 1.5 cm)) inside the nostril and firmly sample the nasal wall at least 4 times. Take approximately 15 seconds to collect the sample. Be sure to collect any nasal drainage that may be present on the swab. Sample both nostril with same swab. And slowly remove swab while rotating it.
- 3) For a nasopharyngeal swab, insert a nasopharyngeal swab into the nostril, swab over the surface of the posterior nasopharynx. Swab should reach depth equal to distance from nostrils to outer opening of the ear. Gently rub and roll the swab. Leave swab in place for several seconds to adsorb secretions. And slowly remove swab while rotating it.


- 4) Place the sample collecting swab into the Extraction Buffer tube containing 300 µL extraction buffer and rotate it more than 5 times to allow extraction.


- 5) Take the sample collecting swab out by pressing and squeezing the sides of the tube to extract the remaining liquid from the swab. Used swab is classified as infectious waste and dispose of used swab properly in accordance with the relevant regulations.

- 6) Press the Dropping Cap onto the Extraction Buffer tube containing the processed sample.

(2) Swab in Viral Transport Media (VTM)


- 1) Mix the specimen in VTM by vortexing.
- 2) Using micro pipette, transfer 300 µL of specimen in VTM to the Extraction bufferTube.
- 3) Press the Dropping Cap onto the Extraction Buffer tube containing the processed sample.

TEST PROCEDURE


1. Open the pouch and take out the Test Cassette. Place it on a flat, dry and clean surface.


2. Invert the Extraction Buffer tube and add 3 drops of processed sample into the sample well on the Test Cassette.

Attention: After the drop is completely absorbed, add the next drop.

3. Read the results in 15~30 minutes after dispensing the sample. Some positive results may appear faster right after the reaction. The result after 30 minutes is invalid.

INTERPRETATION OF RESULTS

- 1. Positive :** If Test line (T) and Control line (C) are appeared in the result window, positive for SARS-CoV-2 antigen.

- 2. Negative :** If only Control line (C) appears in the result window, negative for SARS-CoV-2 antigen.

- 3. Invalid :** If Control line (C) fails to appear, the result is invalid and retest with a new Test Cassette.

QUALITY CONTROL

1. A procedural control is included in the test. A colored line appearing in the control line region (C) is considered an internal procedural control. It confirms sufficient specimen volume, adequate membrane wicking and correct procedural technique.

2. Quality Control materials (positive control swab and negative control swab) can be purchased separately.

LIMITATIONS OF THE SYSTEM

1. The test is for qualitative detection of SARS-CoV-2 antigen in human nasopharyngeal and nasal swab and it does not indicate the quantification of the virus.

2. The test is for *in vitro* diagnostic use only.

3. Negative results do not rule out SARS-CoV-2 infection, particularly in those who have been in contact with the virus. Follow-up testing with a molecular diagnostic should be considered to rule out infection in these individuals.

4. SARS-CoV may cause positive results. SARS-CoV can be detected as a cross reaction.

PERFORMANCE CHARACTERISTICS
1. Limit of Detection : (LoD):

The sensitivity using direct swab is 3.5×10^3 TCID₅₀/mL based on Gamma Irradiated SARS-CoV-2 (BEI Resources, NR-52287, USA-WA1/2020).

The sensitivity using swab in viral transport media is 2.8×10^3 TCID₅₀/mL based on Gamma Irradiated SARS-CoV-2 (BEI Resources, NR-52287, USA-WA1/2020).

2. Cross-Reactivity

SGTi-flex COVID-19 Ag was evaluated with 20 other virus and 12 bacteria. The results show that the SGTi-flex COVID-19 Ag has no cross-reactivity with samples containing tested viruses and bacteria except on SARS-CoV. The results showed no microbial interference with the organisms at the concentrations tested.

Table 1. Microorganisms tested for cross reactivity and interference

No	Virus	Bacteria
1	Alpha Coronavirus (229E)	<i>Hemophilus influenzae</i>
2	Beta Coronavirus OC43	<i>Streptococcus Pneumoniae</i> antigen
3	Human Coronavirus NL63	<i>Candida albicans</i>
4	Beta Coronavirus (MERS) NP protein	Pooled human nasal fluid
5	Beta Coronavirus (SARS-CoV) NP protein	<i>Bordetella pertussis</i>
6	Adenovirus type 5	<i>Mycoplasma pneumoniae</i>
7	Human Metapneumovirus	<i>Chlamydophila pneumoniae</i>
8	Parainfluenza Virus serotype 1	<i>Legionella pneumophila</i>
9	Parainfluenza Virus serotype 2	<i>Staphylococcus epidermidis</i> culture
10	Parainfluenza Virus serotype 3	<i>Mycobacterium tuberculosis</i>
11	Parainfluenza Virus serotype 4	<i>Pneumocystis jirovecii</i> (PJP)
12	Influenza A/H1N1	<i>Staphylococcus aureus</i>
13	Influenza A/H3N2	
14	Influenza A/H5N1	
15	Influenza B	
16	Respiratory Syncytial virus type A	
17	Respiratory Syncytial virus type B	
18	Rhinovirus group A	
19	Human Corona virus HKU1	
20	Enterovirus	

Pneumocystis jirovecii and Human coronavirus HKU1, which were not available for wet testing, were analyzed *in silico* via Basic Local Alignment Search Tool managed by National Center for Biotechnology Information to defined cross-reactivity or interference.

For *Pneumocystis jirovecii*, 45.4% homology was found only a particular part of sequence across 9% of the sequence. It is very unlikely that cross-reaction will occur.

The protein sequences between nucleocapsid protein sequence of human coronavirus HKU1 and nucleocapsid protein of SARS CoV-2 has only 36.7% homology across 82% of sequences. The result of homology between two viruses is relatively very low, but cross-reaction can occur.

3. Analytical Specificity - Interference test

Various concentrations of potential interfering substances were prepared in negative and positive sample. The results show that the SGTi-flex COVID-19 Ag has no interference by the potential interfering substances below which may exist in specimen, such as drugs, and chemical and biological analytes.

Table 2. Interfering substances

No.	Interfering substance	Concentration	No.	Interfering substance	Concentration
1	Albumin	50 mg/ml	15	Tamiflu (Oseltamivir)	6 mg/ml
2	Glucose	1.2 mg/ml	16	Acetaminophen	10 mg/ml
3	Hemoglobin	4 mg/ml	17	Ibuprofen	5 mg/ml
4	Bilirubin	5 mg/ml	18	Aspirin	2 mg/mL
5	mucin	1.0 %	19	Naso GEL	5% v/v
6	Whole blood	4.0 %	20	Oxymetazoline	0.1 mg/mL
7	Phenylephrine hydrochloride	10 mg/ml	21	Cromolyn	0.03 mg/mL
8	Dexamethasone	0.6 mg/ml	22	Zicam	5% v/v
9	Flunisolide	2.5 mg/ml	23	Alkalol	10% v/v
10	Budesonide	1 mg/ml	24	Mupirocin	10 mg/mL
11	Benzocaine	5 mg/ml	25	Fluticasone Propionate	5% v/v
12	Menthol	40 mg/ml	26	Sore Throat Phenol Spray	15% v/v
13	Zanamivir	10 mg/ml	27	Heparin sodium salt	3000 U/L
14	Tobramycin	20 mg/ml			

4. High-dose Hook Effect

There is no hook effect was observed at high levels of Gamma Irradiated SARS-CoV-2 (BEI Resources,

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