SGTi-flex COVID-19 Ag

IVD REF CAGT025E0, CAGT025E1 CE

English

INTENDED USE

SGTi-flex COVID-19 Ag is an immunoassay for qualitative detection of SARS-CoV-2 antigens directly from nasopharyngeal and oropharyngeal swab specimens. The test is used as an aid in the rapid diagnosis of SARS-CoV-2 viral infections.

SUMMARY AND EXPLANATION

The novel coronavirus (SARS-CoV-2) was identified in December 2019, and in February 2020, the World Health Organization (WHO) officially named the disease caused by SARS-CoV-2 as COVID-19 (Coronavirus Disease 2019). Belonging to the family Coronaviridae, it has a positive-sense single-stranded RNA and can be transmitted between people. The coronaviruses identified for human infection include 229E, NL63 belonging to α-Coronaviruses and HKU1, OC43, SARS-CoV, MERS-CoV belonging to β-Coronaviruses.

The new coronavirus was published under the name of SARS-CoV-2, with 80% of genetic similarity to SARS-CoV by ICTV International Committee on Taxonomy of Viruses.

COVID-19 spreads mainly through respiratory droplets, which cause lethargy, fever, dry cough, and dyspnea when infected. It can be even led to death with its severe symptoms like sepsis, MOF (Multiple Organ Failure) and ARDS (Acute Respiratory Distress Syndrome). It is more contagious than SARS which caused more than 800 deaths and 8,000 infected patients. Moreover, it has an incubation period of about 3 days to up to 16 days and becomes a big threat as infectivity appears even during the incubation period. There is currently no specific treatment for COVID-19, and rapid and accurate diagnosis is an important issue for isolation of patients with symptoms of suspected COVID-19.

PRINCIPLE

SGTi-flex COVID-19 Ag is an immunoassay for qualitative detection of SARS-CoV-2 antigens directly from nasopharyngeal and oropharyngeal swab specimens. The SARS-CoV-2 antigens are extracted from swab in the extraction buffer and the extracted sample solutions are loaded to the sample well of the Test Cassette. When the sample is loaded, the detection antibody binds to SARS-CoV-2 antigen and flows through the membrane. The detection antibody-gold conjugate and SARS-CoV-2 antigen move to the test line area and are accumulated by the capture antibody immobilized on the membrane. This leads to the formation of a reddish colored band. The intensity of the band depends on quantity of SARS-CoV-2 antigen and the test results are interpreted by user's eye according to the instructions for use.

MATERIALS SUPPLIED

SGTi-flex COVID-19 Ag has two types (CAGT025E0, CAGT025E1) that differ in components. Please see the REF number described on the label of the packaging.

CAGT025E0		CAGT025E1	
Test Cassette	25	Test Cassette	25
Sample Collection swab	25	Sample Collection swab	25
Dropping Cap	25	Dropping Cap	25
Extraction Buffer	25 (0.3 mL/Tube)	Extraction Buffer	1 (10 mL/Bottle)
Instructions for Use	1	Instructions for Use	1

STORAGE AND STABILITY

• Store SGTi-flex COVID-19 Ag Test Cassette and Extraction Buffer at 2–30°C (36–86°F).
• If SGTi-flex COVID-19 Ag Test Cassette and Extraction Buffer are stored in cold storage, allow them for 30 minutes to return to room temperature before testing.
• Do not open the pouch of Test Cassette until ready to use. After opening aluminum pouch, Test Cassette should be used immediately.
• Keep away from direct sunlight.

WARNING AND PRECAUTIONS

- For in-vitro diagnostic use only.
- For use by trained laboratory personnel or healthcare professionals. The result of this test should not be the sole basis for the diagnosis. Confirmatory testing is required.
- Clinical diagnosis through this product should be made through a comprehensive review of the specialist based on other test methods and clinical symptoms.
- Please read the instruction carefully before you begin the test and follow the procedure correctly.
- It is prohibited to reuse Test Cassettes because they are single use only.
- The test result after the expiry date is not reliable.
- Test Cassette is sensitive to moisture and should be stored in a sealed pouch until use. Use Test Cassette immediately after opening the pouch.
- Do not use the Test Cassette if it is broken or the pouch is not stored in sealed.
- Samples and Test Cassette must be at room temperature before testing.
- It is an in-vitro diagnostic product and the risk of infection is low because there is no direct contact with the human body. However please be cautious when handling Test Cassette and samples because of the use of clinical samples containing potential infectious sources. Dispose of the used samples and Test Cassettes properly in accordance with the relevant regulations.
- Smoking and eating are prohibited at test site when handing specimens or kit reagents.

TEST PREPARATION

1. **Test should be done immediately after sample collecting.**
 - 1) If sample swabs are not used immediately after sample collection, specimen is recommended to be stored in deep freezer at -70°C (or in dry ice or liquid nitrogen). A freezer at -20°C is NOT recommended, if the specimen is stored at 2-8°C, it can be stored up to 72 hours.
2. **Preparation before Test**
 - 1) All samples and reagents should be stored at room temperature and stayed homogenous 15–30 minutes prior to testing.
 - 2) Test cassette is moisture sensitive so should be used **immediately** after opening.

SAMPLE COLLECTION

SGTi-flex COVID-19 Ag can be performed with nasopharyngeal swab and oropharyngeal swab.


1. Before sample collection, prepare the extraction buffer.
- 1) For CAGT025E1, dispense Extraction Buffer into the Extraction Tube until it flows up to the fill-line (300 µL).
- 2) For CAGT025E0, remove the sealing foil from the Extraction Buffer tube.
2. Place the Extraction Tube in the tube rack.
3. SGTi-flex COVID-19 Ag uses the sample of nasopharyngeal swab and oropharyngeal swab.
- 1) Please use single use sample collecting swab.
- 2) Insert a nasopharyngeal swab into the nostril of the patient, swab over the surface of the posterior nasopharynx. Swab should reach depth equal to distance from nostrils to outer opening of the ear.

※ The sample collection swab provided by SGTi-flex COVID-19 Ag is used for nasopharyngeal swab.

3) Gently rub and roll the swab. Leave swab in place for several seconds to absorb secretions. And slowly remove swab while rotating it.

4) Place the sample collecting swab into the Extraction Buffer tube containing 300 µL extraction buffer and rotate it more than 5 times to allow extraction.

5) Take the sample collecting swab out by pressing and squeezing the sides of the tube to extract the remaining liquid from the swab. Used swab is classified as infectious waste and dispose of used swab properly in accordance with the relevant regulations.

<nasopharyngeal swab>

x5

20~30min

- Various concentrations of potential interfering substances were prepared in negative and positive sample. The results show that the SGTi-flex COVID-19 Ag has no interference by the potential interfering substances below which may exist in specimen, such as prescription/OTC drugs, and elevated levels of chemical and biological analytes.

Table 3. Interfering substances

No. Interfering substance Concentration No. Interfering substance Concentration

1 Albumin 50 mg/mL 10 Menthol 40 mg/mL

2 Glucose 1.2 mg/mL 11 Zanamivir 10 mg/mL

3 Phenylephrine hydrochloride 10 mg/mL 12 mucin 1.0 %

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LIMITACIONES DEL SISTEMA

1. The test is for qualitative detection of SARS-CoV-2 antigen in human nasopharyngeal and oropharyngeal swab and it does not indicate the quantification of the virus.
2. The test is for in-vitro diagnostic use only.
3. Negative results do not rule out SARS-CoV-2 infection, particularly in those who have been in contact with the virus. Follow-up testing with a molecular diagnostic should be considered to rule out infection in these individuals.
4. SARS-CoV may cause positive results. SARS-CoV can be detected as a cross reaction.

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